

Highlights of the two-piece Z5c implant

Z5c-40-10


> Implant 4 mm × 10 mm

Same outer contour and surface like monotype Z-SYSTEMS implants ensures transferable Results:

- same drilling set
- excellent osseointegration with SLM® laser surface
- very good long-term results

Experience from over 60,000 implants placed

- Self-tapping thread with bone chip reservoir reduces insertion torque
- Excellent osseointegration through proven, patented proprietary SLM® laser surface technology and gentle plasma sterilization



Can be prepped directly in situ: Abutment and implant shoulder

- individual aesthetics
- Fast and cost-effective restorations – like a natural tooth

- Perfect individual aesthetics thanks to preppable preparation margin.

- The implant connects a conical basic shape with a dynamic self-tapping thread.



Z-SYSTEMS
ceramic implants